

# Skin Closures

Butterfly Closures

Indermil

Steri-Strips

Sutures

Davis & Geck

Ethicon

US Surgical

Skin Staplers

Staple Removers/Extractors

Suture Removal Trays

Non-Absorbable Sutures

**NYLON** – Synthetic uncoated monofilament. Smooth tissue passage, minimal tissue reaction. General skin closures, plastic surgery, and removable subcuticular.

AUGHT SIZE	ETHICON - ETHILON			DAVIS & GECK - DERMALON			US SURGICAL - MONOSOF		
	NEEDLE TYPE	PRODUCT CODE	BOX QTY	NEEDLE TYPE	PRODUCT COCE	BOX QTY	NEEDLE TYPE	PRODUCT CODE	BOX QTY
3 / 0	FS-1	130-P663	12	CE-6	111-1644-43	12	C-14	111-SN-663G	12
3 / 0				CE-6	111-1644-41	36	C-14	111-SN-663	36
3 / 0				CE-4	111-1656-41	36	C-13	111-SN-653	36
4 / 0	FS-2	130-P662	12	CE-4	111-1656-33	12	C-14	111-SN-662G	12
4 / 0				CE-4	111-1656-31	36	C-14	111-SN-662	36
4 / 0				CE-2	111-1641-31	36	C-1	111-SN-872	36
4 / 0	P-3	130-P699	12	PRE-2	111-1807-33	12	P-13	111-SN-1699G	12
5 / 0	FS-2	130-P661	12	CE-4	111-1656-23	12	C-13	111-SN-661G	12
5 / 0				CE-4	111-1656-21	36	C-13	111-SN-661	36
5 / 0	C-2	130-P668	12	CE-2	111-1641-21	36	C-1	111-SN-871	36
5 / 0	P-3	130-P698	12	PRE-2	111-1807-23	12	P-13	111-SN-3698	36
6 / 0	C-2	130-P667	12	CE-2	111-1641-11	36	C-1	111-SN-667	36
6 / 0				CE-4	111-1656-11	36	C-13	111-SN-650	36
6 / 0	FS-3	130-P660	12	CE-3	111-1647-11	36	C-12	111-SN-660	36
6 / 0	P-1	130-P697	12	PRE-1	111-1804-13	12	P-10	111-SN-3697	36
6 / 0	P-3	130-1698G	12	PRE-2	111-1807-13	12	P-13	111-SN-3696	36

**POLYPROPYLENE** – Linear hydrocarbon polymer monofilament. Excellent tensile strength. Smooth with minimal tissue reactivity. General closure, skin, vascular, cardiovascular orthopedics, pediatric surgery, ophthalmology, plastic and microvascular surgery.

AUGHT SIZE	ETHICON - PROLENE			DAVIS & GECK - SURGILENE			US SURGICAL - SURGIPRO		
	NEEDLE TYPE	PRODUCT CODE	BOX QTY	NEEDLE TYPE	PRODUCT COCE	BOX QTY	NEEDLE TYPE	PRODUCT CODE	BOX QTY
2 / 0	FS	130-8685H	36						
3 / 0	FS-2	130-8665H	36	CE-4	111-3420-43	12	C-13	111-SP-665G	12
4 / 0	FS-2	130-8683G	12	CE-4	111-3420-33	12	C-13	111-SP-683G	12
5 / 0	FS-2	130-8661G	12	CE-4	111-3420-23	12	C-13	111-SP-661G	12
6 / 0	P-1	130-8697G	12	PRE-1	111-3410-13	12	P-10	111-SP-1697G	12

**SILK** – Braided, siliconized proteinaceous thread spun of silkworm larval. Soft, supple and smooth. Good knot security. Not as strong as comparable sizes of synthetic materials. Excellent handling characteristics. Used on G.I. tract, neurosurgery, general closure, skin, plastic, ophthalmic and cardiovascular surgery.

AUGHT SIZE	ETHICON – SILK			DAVIS & GECK – SILK			US SURGICAL - SOFSILK		
	NEEDLE TYPE	PRODUCT CODE	BOX QTY	NEEDLE TYPE	PRODUCT COCE	BOX QTY	NEEDLE TYPE	PRODUCT CODE	BOX QTY
2 / 0	FS	130-685G	12						
3 / 0	FS-1	130-684G	12				C-14	111-SS-684G	12
4 / 0				CE-4	111-1074-41	36	C-13	111-SS-694	36
5 / 0	FS-2	130-682G	12	CE-4	111-1074-23	12	C-13	111-SS-682	36

Absorbable Sutures

**GUT** – Twisted proteinaceous sheep or cow intestinal lining. Good tensile strength. Absorption is by phagocytosis, which is unpredictable, with gradual loss of strength. Suture is wiry if allowed to dry. Gut can be chromacized or untreated. General closure, OB/GYN, ophthalmology, urology, orthopedics, bowel anastomosis, and ligation. Not used where long-term tissue approximation is required. Faster absorption in infected tissue.

	ETHICON – PLAIN GUT			DAVIS & GECK – PLAIN GUT			US SURGICAL - SURGIGUT		
AUGHT SIZE	NEEDLE TYPE	PRODUCT CODE	BOX QTY	NEEDLE TYPE	PRODUCT COCE	BOX QTY	NEEDLE TYPE	PRODUCT CODE	BOX QTY
3 / 0	FS-1	130-H823H	36	CE-6	111-68-51	36	C-14	111-SG-823	36
3 / 0	FS-2	130-H822G	12	CE-4	111-65-41	36	CE-4	111-SG-822	12
4 / 0	FS-2	130-H821G	12	CE-4	111-65-31	36	CE-4	111-SG-752	36
5 / 0	FS-3	130-686G	12	PRE-2	111-74-23	12	PRE-2	111-SG-3696	12

	ETHICON – CHROMIC GUT			DAVIS & GECK – CHROMIC GUT			US SURGICAL – CHROMIC GUT		
AUGHT SIZE	NEEDLE TYPE	PRODUCT CODE	BOX QTY	NEEDLE TYPE	PRODUCT COCE	BOX QTY	NEEDLE TYPE	PRODUCT CODE	BOX QTY
3 / 0	FS-2	130-636G	12	CE-4	111-566-41	36	C-13	111-CG-676	36
4 / 0	FS-2	130-635G	12	CE-4	111-566-31	36	C-13	111-SG-675	36


**POLYGLYCONATE** – Synthetic monofilament, non-antigenic and non-pyrogenic. Smooth tissue passage. Minimized need for knot repositioning. For use in all types of soft tissue approximation, G.I. anastomosis, facial closures and caesarean section. Not for use in cardiovascular, ophthalmic or microsurgery.

	ETHICON – PDS II			DAVIS & GECK – DEXON II			US SURGICAL		
AUGHT SIZE	NEEDLE TYPE	PRODUCT CODE	BOX QTY	NEEDLE TYPE	PRODUCT COCE	BOX QTY	NEEDLE TYPE	PRODUCT CODE	BOX QTY
3 / 0				CE-6	111-6618-41	36			
3 / 0	SH	130-7316H	36	T-5	111-6233-41	36			
4 / 0	FS-2	130-Y422H	36	CE-4	111-6608-31	36			
5 / 0				CE-4	111-6608-21	36			

**POLYGLYCOLIC ACID** – Braided, polycaprolate non-toxic coating for improved handling and tissue passage. Used for most procedures where absorption is required. Not to be used where long-term tissue approximation is required

	ETHICON – VICRYL			DAVIS & GECK – MAXON			US SURGICAL		
AUGHT SIZE	NEEDLE TYPE	PRODUCT CODE	BOX QTY	NEEDLE TYPE	PRODUCT COCE	BOX QTY	NEEDLE TYPE	PRODUCT CODE	BOX QTY
2 / 0	SH	130-J317H	36						
3 / 0	FS-1	130-J442H	36	CE-6	111-9618-41	36			
3 / 0	FS-2	130-J423H	36	CE-4	111-9608-41	36			
4 / 0	FS-2	130-J422H	36	CE-4	111-9608-31	36			
5 / 0	FS-2	130-J421H	36	CE-4	111-9608-21	36			
5 / 0	P-3	130-J493G	12	PRE-2	111-9536-*21	36			
6 / 0	P-3	130-J492G	12	PRE-2	111-9536-11	36			

Skin Closures



Skin Adhesive  
Translucent Medical Adhesive For Closing  
Wounds and Lacerations.

Indermil, 5 x 0.5gm, 12/bx

111.IND012



3M Steri-Strips  
Reinforced Skin Closures.

1/8 x 3", 5 x 50/bx, 4bx/cs

105.R-1540

1/4" x 1 1/2", 5 x 50/bx, 4bx/cs

105.R-1542

1/4 x 3", 3 x 50/bx, 4bx/cs

105.R-1541

1/4 x 4", 10 x 50/bx, 4bx/cs

105.R-1546

1/2 x 4", 6 x 50/bx, 4bx/cs

105.R-1547

Steri-Strip Organizer / Dispenser, ea

600.MMM-1555



Butterfly Wound Closures

Medium, 100/bx

501.1975

Large, 100/bx

501.1990


Suture / Staple Removal



Skin Staplers  
Sterile, Disposable Single Use.

3M, 15 Shot, 12/bx

105.DS-15



Staple Remover / Extractor  
Sterile, Disposable Single Use.

Proximate, 12/bx

130.PSX

Appose, 12/bx

111.8039-12

3M Tweezer Style, 30/cs

105.SR-1

3M Scissor Style, 30/cs

105.SR-3

**\*For Reusable Staple Removers See  
Section 14 – Instruments.**



Suture Removal Trays

Contents: 1- 3 x 3 Gauze Sponge,  
1- Metal Forceps, 1- Suture Scissors,  
1- CSR Wrap. ea, cs/50

600.MAI-69239