



RESPIRATORY




| | |
|-------------------------------------|---|
| Catheter Introducing Forceps | 2 |
| CPR Barrier | 3 |
| CPR Masks | 3 |
| Endotracheal Tubes | 2 |
| Guedel Airways | 3 |
| Laryngeal Airway Masks | 2 |
| Laryngoscope Kits | 3 |
| Manual Pulmonary Resuscitator | 3 |
| Oropharyngeal Airways | 2 |
| Non-Rebreather Masks | 2 |
| Oxygen Cannulas | 1 |
| Oxygen Masks | 1 |
| Oxygen Tank Cart | 1 |
| Oxygen Tubing | 1 |



OXYGEN CANNULAS, OXYGEN MASKS & NON-REBREATHER MASKS

Oxygen Cannulas




| Description | Size | Quantity | Order Number |
|--|-------|----------|-----------------|
| Amsino AMsure Oxygen Cannula Flared prongs, over the ear design, with 7 ft kink-resistant tubing | Adult | Each | 600.AMS-AS75085 |

Oxygen Mask



| Description | Size | Quantity | Order Number |
|---|-------|----------|-----------------|
| Amsino AMsure , with 7 foot tubing | Adult | Each | 600.AMS-AS74010 |

Oxygen Tubing



| Description | Size | Quantity | Order Number |
|-----------------------------------|------|----------|-----------------|
| Amsino Oxygen Tubing Only, 7 foot | | Each | 600.AMS-AS76007 |

Non-Rebreather Oxygen Mask



| Description | Size | Quantity | Order Number |
|---|-----------|----------|-----------------|
| Amsino Non-rebreather Oxygen Mask, (2) one way valves, with bag and 7 foot tubing | Adult | Each | 600.AMS-AS75010 |
| Partial Non-rebreather Oxygen Mask, (1) one way valve with 7 foot tubing | Pediatric | Each | 103.1098 |

A non-rebreather mask is a device used in medicine to assist in the delivery of oxygen therapy. A non-rebreather mask requires that the patient can breathe unassisted, but unlike low-flow nasal cannulae, the non-re-breather mask allows for the delivery of higher concentrations of oxygen. An ideal non-rebreather mask does not permit air from the surrounding environment to be inhaled. This type of mask also has a reservoir bag attached



Oxygen Tank Cart

| Description | Quantity | Order Number |
|--|----------|--------------|
| Dynarex Single Cylinder Oxygen Tank Cart | Each | 609.10556 |

AIRWAY MANAGEMENT



Catheter Introducing Forceps

| Description | Size | Quantity | Order Number |
|--|-------|----------|--------------|
| Magill Catheter Introducing Forceps, 9-3/4", Magna | Adult | Each | 115.M40-0120 |



Laryngeal Airway masks

| Description | Size | Quantity | Order Number |
|--|----------------|----------|-----------------|
| Smiths-Medical Portex Soft seal | Infant/Child | 2 Each | 103.100/222/200 |
| | Child | 3 Each | 103.100/222/300 |
| | Adult 50-70kg | 4 Each | 103.100/222/400 |
| | Adult 70-100kg | 5 Each | 103.100/222/500 |

Clear tubing. Atraumatic tip. Non-return valve. Printed pilot balloon. 15mm connector. Silicone. Latex-Free. Disposable.



Oropharyngeal Airways

| Description | Size | Quantity | Order Number |
|-------------|------|----------|--------------|
| Rusch | | Set/7 | 103.123325 |

Soft PVC construction • Inner striations for smooth suction catheter insertion • Beveled tip for comfortable,



Endotracheal Tube

| Description | Size | Quantity | Order Number |
|---|-------|----------|--------------|
| Hi-lo Tracheal tube, cuffed Murphy hooked eye tip | 6.0mm | Each | 222.TY86109 |
| | 6.5mm | Each | 222.TY86110 |
| | 7.0mm | Each | 222.TY86111 |
| | 7.5mm | Each | 222.TY86112 |
| | 8.0mm | Each | 222.TY86113 |

Pilot balloon and mechanical self-sealing valve. Large cuff diameter of 1-1/2 times the average trachea maintains proper low-pressure seal during tracheal expansion and contraction. Hi-Lo® offers better protection against aspiration with its large cuff resting volume and long length to reduce pressure against tracheal tissue.

RESPIRATORY



AIRWAY MANAGEMENT



Laryngoscope Kit

| Description | Size | Quantity | Order Number |
|--|------|----------|--------------|
| ADC Satin Laryngoscope Set | | | 617.4079 |
| Includes: #1,2,3 and 4 Macintosh Blades. Medium C handle, Small AA Handle, Fitted hard case. | | | |

Guedel Airways (Oropharyngeal)

| Description | Quantity | Order Number |
|------------------------------------|----------|--------------|
| Contact us for further information | 8/Set | 5000 |

RESUSCITATION



Manual Pulmonary Resuscitator (BVM)

| Description | Size | Quantity | Order Number |
|------------------------------|-----------|----------|--------------|
| 2500cc/ml bag, 10 ft tubing, | Pediatric | Each | 609.5804 |
| | Adult | Each | 609.5808 |

A large rubber-like bag with oxygen tubing, duck-bill valve and exhalation port, reservoir tube or reservoir bag and other components.

CPR Resuscitator Mask

| Description | Quantity | Order Number |
|---|----------|--------------|
| ADC CPR valve mask resuscitator in case | Each | 600.ADC-4053 |

High efficiency and low resistance one-way valve minimizes the possibility of cross-contamination. 15mm O.D. connector can be attached to any standard ventilation and oxygen equipment. For use on adults, and inverted for infants/children. Ergonomically shaped mask with soft-air cushion for leak-free performance, easy grasp, and patient comfort Disposable valve to avoid cross contamination Elastic green strap to keep mask in place on patient. Not made with natural rubber latex. Operating instructions, disposable gloves not made with natural rubber latex, and alcohol wipe included. Pocket storage case.



CPR Barrier

| Description | Quantity | Order Number |
|-----------------------------------|----------|--------------|
| Adsafe CPR Face Shield | Each | 617.4055W |
| Adsafe CPR Face Shield with valve | Each | 617.4056W |

User friendly design protects the rescuer and helps eliminate hesitation when performing CPR. Semi-transparent face shield allows visualization of patient's lips, mouth, and nose while protecting rescuer. Printed instructions and orientation guides right on the shield help ensure correct use, even by those with little or no CPR experience Filter helps protect the rescuer. Single use, Not made with natural rubber latex. Measures 12"L x 8"W when open, a miniscule 1 7/8" diameter when stored in its case. Comes in compact hard cylindrical case with handy key ring clasp



RESPIRATORY